eClinicalWorks

eCW Compatibility with Google Chrome 80

With the next major version of Google[™] Chrome[™], Chrome 80, Google is changing the handling of browser cookies. This change will prevent certain eClinicalWorks[®] application windows from functioning as expected.


Google plans to release Chrome 80 starting Tuesday, February 4, 2020. eClinicalWorks is working to address these concerns but cannot guarantee a seamless solution will be in place before Google distributes the new release to client machines.

Required actions:

Chrome 80 is not supported for use with certain component products at this time, and eClinicalWorks recommends that end users do not upgrade to Chrome 80.

Clients that have already upgraded to Chrome 80 can configure a setting to restore browser functionality:

- 1. Enter the following in the Chrome address bar: chrome://flags/#same-site-by-default-cookies
- 2. Next to the SameSite by default cookies setting, select Disabled from the drop-down list:



Questions, issues, or requests:

Create a support case on the my.eclinicalworks.com Customer Portal with any questions regarding the operation or deployment of this release.

APPENDIX A: DOCUMENTATION TERMS AND CONDITIONS

Documentation: This document contains information that is confidential and proprietary to eClinicalWorks, LLC and is intended for use solely by its authorized licensees ("Authorized Entities"). In accordance with these conditions and contractual agreements by the user, this document may not be copied, displayed, distributed, published, or otherwise reproduced, otherwise used, transmitted, or in any form or otherwise made available or used by anyone other than the authorized client to whom this document was originally delivered without the prior written consent of eClinicalWorks, LLC. Pursuant to its agreement, Authorized Entities may receive copies of or access to certain written technical support or explanatory documents regarding eClinicalWorks' software; eClinicalWorks' services; and/or eClinicalWorks' internal policies and procedures (collectively, "Documentation").

The examples, images, and scenarios presented in this documentation are solely for explanatory use of the software and its functionality and should not be construed as directives for clinical or medical decisions; the user is ultimately and completely responsible for clinical and/or medical decisions made regarding patient care.

eClinicalWorks documentation may contain hyperlinks to external sites and/or third-party vendors for functional, informational, or instructional purposes. Use of these external links is at the user's and licensee's risk.

eClinicalWorks assumes no responsibility for errors or omissions that may appear in this publication and reserves the right to change this publication at any time without notice. *All users are instructed to consult the latest version of the Documentation at my.eclinicalworks.com for the latest updates to the Documentation at all times.*

Once obtained, the distribution or posting this proprietary document on the Internet for public and/or private use is strictly prohibited. This restriction includes Internet websites, forums, blogs, private or public portals, or any other electronic means of sharing beyond the intended, licensed user.

APPENDIX B: NOTICES

The following appendix lists the trademark and copyright information for this document.

Trademarks

eClinicalWorks[®]

eClinicalWorks[®] is a registered trademark of eClinicalWorks, LLC.

All other trademarks or service marks contained herein are the property of their respective owners.

Google™

Chrome™

Google and Chrome are trademarks, registered trademarks, or service marks of Google, Inc.